Profin® PDH

Multilingual A: EN, SV, DE, ES, FR, IT, BG, CS, DK, EL

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BL-6734-00_2017-04-05

EN - Profin PDH Instruction for use – maintenance information

Prior to initial use - always read the instructions for use. Only for USA: prescription use – for use by dental clinicians only.

Indication for use Finish and polish tooth surfaces.

Precaution Insert the shank of the tip into the holder. Check that the tip fits safely by pulling the tip gently with your fingers. A firm resistance should be sensed.

Caution If a tip comes out too easily from the holder, discard and use another tip.

Technical description	Recommendation of use
Holder	
Profin PDH angled end	For use in distal areas
Profin PDH straight end	For use in anterior areas
Variation of mode	
A. Tip in a self-steering mode (The tip can rotate along its axis.)	For following the contour of the tooth easily. Normally used during polishing and finishing.
B. Tip in fixed mode in one of the indexed positions. Eight positions/settings.	For improved tactility/accuracy. E.g. for creating individual structures for aesthetic appearance. - To change the position – pull out the tip slightly, turn into desired position and push it back into a suitable notch in the holder.
Durability of tips	
Diamond or Tungsten Carbide coated	The durability depends on the circumstances, for continued efficiency it is recommended not to use them more than 30-40 times.

Recommendations

- A general recommendation is to use a more coarse tip when larger amount of material is to be removed. Removal of material should always be made successively from tips with more coarse grit to finer.
- Do not change directly from super-coarse to super-fine tips in one step, but use intermediate m coarse-m edium.
- The black Tungsten-Carbide-coated tips can always be used for finishing when there is a risk for scarring enamel, since Tungsten-Carbide does not scar or scuff enamel.

Maintenance, Profin PDH

Limitations of reprocessing	Check for wear of holders
1. Preparation at the point of use	None.
2. Decontamination	
Cleaning	Manual Flush off larger debris under running water. Use a not-metal brush to avoid scratching the instrument. Wipe off with a soft cloth. Without or with mild detergents. Always follow the manufacturer's instruction.
Disinfection	Without or with mild detergents. Always follow the manufacturer's instruction.
Drying	Maximum 95 °C (203 °F). Without or with mild detergents. Always follow the manufacturer's instruction.
Inspection	Without or with mild detergents. Always follow the manufacturer's instruction.
3. Sterilization - optional	
	134°C Place in sterilization bag of paper/clear film (according to applicable standards). Check correct sealing of the bag, 134 °C (273 °F)
Storage	Protect from dust and dirt.

LTA-tips

Limitations of reprocessing	See Durability above.
1. Preparation at the point of use	Check for wear of tip
2. Decontamination	
Cleaning	Without or with mild detergents. Always follow the manufacturer's instruction.
Disinfection	Manual Rub against the cleaning stone to remove clogged debris, flush under running water.
Drying	Thermal disinfectors Maximum 95 °C (203 °F). Without or with mild detergents. Always follow the manufacturer's instruction.
Inspection	Let air dry. Visually inspect for appearance and function. Eliminate worn tips.
3. Sterilization - optional	
	Place in sterilization bag of paper/clear film (according to applicable standards). Check correct sealing of the bag, 134 °C (273 °F)
Storage	Protect from dust and dirt.

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