

AG-3 accessories



AG-3 DA

upper and lower model jaws combined in a demonstration articulator allowing simple functional movements



AG-3 WOK

elastic removable gingiva for upper jaw, positioned over the model base, held in place by the model teeth suitable for AG-3 OK



AG-3 WUK

elastic removable gingiva for lower jaw, positioned over the model base, held in place by the model teeth suitable for AG-3 UK



AG-3 WOK 20

elastic removable gingiva for upper jaw, similar to AG-3 WOK, however the gingival contour is taller to represent that of a 20-year-old patient

AG-3 WUK 20

elastic removable gingiva for lower jaw, similar to AG-3 WUK, however the gingival contour is taller to represent that of a 20-year-old patient



AG-3 WOKH

By changing the gingivae, AG-3 models can be transformed into AG-3 H dental hygiene models as shown on page 51. Further accessories can be found there.

AG-3 WUKH

elastic removable gingiva for oral hygiene model upper jaw, held in place by the model teeth

elastic removable gingiva for oral hygiene model lower jaw, held in place by the model teeth