

MY LOAD IS ...	MY LOAD CONTAINS ...	MY LOAD IS MADE OF ...	THEN I CHOOSE PROGRAM ...		THAT LASTS APPROX. ...
Wrapped or both wrapped and unwrapped	Only narrow lumen instruments or a mix of instruments	Metal only	No. 4	B134 Hollow Wrapped	29 min
		Mixed materials	No. 3	B121 Hollow Wrapped	40 min
	Only porous material	Heat resistant porous material	No. 9	B134 Porous	40 min
		Heat sensitive porous material	No. 8	B121 Porous	50 min
	Only solid and/or simple hollow instruments	Metal only	No. 6	S134 Solids Wrapped	25 min
		Mixed materials	No. 5	S121 Solids Wrapped	37 min
Unwrapped only	Only narrow lumen instruments or a mix of instruments	Metal only	No. 12	B134 Hollow Open	26 min
		Mixed materials	No. 3	B121 Hollow Wrapped	40 min
	Only porous material	Heat resistant porous material	No. 9	B134 Porous	40 min
		Heat sensitive porous material	No. 8	B121 Porous	50 min
	Only solid and/or simple hollow instruments	Metal only	No. 11	S134 Rapid	20 min
		Mixed materials	No. 10	S121 Rapid	31 min
Related to Creutzfeldt-Jakob disease	A mix of instruments	Mixed materials	No. 7	B134 Prion	44 min