

DENTATUS PROFIN PDH

THE ULTIMATE INSTRUMENT FOR MANUAL FINISHING

Safe and precise dental finishing, even in the most inaccessible areas!

The Advantages

- Contra angled for easy posterior and buccal access
- Interproximal adjustments
- Overhang reduction
- Marginal adjustments
- Professional tooth cleaning

Profin PDH is an ergonomic hand instrument of the highest quality based on the established Profin/EVA system. Profin PDH enables you to reach, access and finish areas that are otherwise hard to reach with traditional diamond burs, discs or strips.

The interchangeable Lamineer diamond tips are available in various grits and lengths for safe and precise manual finishing.

The hand instrument is also an optimal choice for easy and comfortable professional mechanical tooth cleaning with its contra angled part. The ergonomic design facilitates the work of the care provider and highly increases patient comfort.

Contra angled for posterior usage

Dual modalities: fixed or self-stering tip


Straight for anterior usage



Dentatus Profin PDH for Different Areas of Use



Straight part for anterior fine finishing of fillings for optimal esthetics



To reach, access and remove interproximal overhangs of cement, composite and amalgam



Contra angled for easy posterior access



Interproximal professional mechanical tooth cleaning



Dentatus AB
www.dentatus.com
SE: +46 8 546 509 00
info@dentatus.se
US: +1 212 481 1010
dentatus@dentatus.com

