

Achieve the highest levels of dental finishing!

Based on the established Profin/EVA system, Dentatus is proud to present the new Profin PDH, a high-quality hand instrument with a wide range of interchangeable diamond tips for safe and precise finishing work.

Profin PDH enabling you to achieve the highest levels of dental finishing during daily use, even in the most inaccessible areas.

IDEAL FOR:


- Interproximal adjustments
- Marginal adjustments
- Overhang reduction



Straight for anterior usage

Dual modalities: fixed/locked or self-steering tip

Contra angled for posterior usage



Profin PDH is perfect for creating optimal shapes and embrasures, giving heightened aesthetic results and increased patient satisfaction. With ergonomic, tactile design and tough anodized, scratch-resistant surface for optimal durability.

DENTATUS LAMINEER TIPS FOR PROFIN/EVA

The Lamineer Tips remove overhangs in all leading dental filling materials, including composite, amalgam, and cement. Access and refinish residual bonding material in inaccessible and sensitive gingival and interproximal areas.

- Eliminate ditching associated with rotating instruments and imprecise strip work
- High quality diamond tips in a broad range of grain sizes
- Fixed or self-steering tip settings and world-renowned reciprocating movement for optimal preparation and access



The spatula and knife-edge Lamineer tips are safe-sided, and are available in a wide range of abrasive surfaces and grits. The disposable plastic tips are designed for interproximal polishing with a your choice of polishing compounds. No hygiene issues: Lamineer tips are fully autoclavable



PROFIN PDH/SK CONTAINS:

- 1 x PDH
- 1 x LTA-150
- 1 x LTA-75

INDISPENSABLE DURING CAD-CAM PREPARATION

Whether you are preparing for an in-lay or on-lay, refining vertical proximal walls and margins or finishing interproximal contours, you will find Profin's precision an indispensable asset to your CAD-CAM restorations.



Dentatus AB
Sweden +46 8 546 509 00, info@dentatus.se

MADE IN SWEDEN

www.dentatus.com
USA + 1 212 481 1010, dentatus@dentatus.com