



## sundries

### frasaco – crowns and more

#### frasaco strip crowns:

- renowned moulds for fast and clean temporary restorations
- elastic, thin and easy to cut, packed antistatically
- pedodontal strip crown available for restoration of badly damaged deciduous teeth
- wide range covering all required sizes
- assortments for proper storage and choice, refill-packages available

#### frasaco patented dispensers:

- quick and simple single-handed operation
- band exactly in the required length with neat and even cuts
- economical use and reduced waste
- one time purchase of the dispenser, cost-saving refill rolls available



Medizinprodukt Klasse 1



frasaco strip crowns,  
matrix strips and ma-  
trix bands are medical  
devices and therefore  
intended for unique  
use only

# matrix strips and dispensers

Tear resistant, plastic universal strips for use with auto-curing and light-curing composites and cements.  
To be used when shaping the restorations for class II, III, IV and V cavities

- Transparent
- Extra thin or standard thickness (retains its stability when being used for larger cavities)

**Note: the part of the transparent matrix strips with the colour mark must not enter the mouth of the patient.**

Enjoy the one time purchase of the patented matrix band dispenser, refill rolls are available at an attractive price

universal strips, straight, length: 120 mm, width: 8.5 mm, package containing 50 pieces

- |                 |                                            |
|-----------------|--------------------------------------------|
| <b>US-120KE</b> | transparent<br>extra thin: 0.05 mm         |
| <b>US-120KN</b> | transparent<br>standard thickness: 0.08 mm |



universal strips, curved, length: 75 mm, width: 8.5 mm, package containing 50 pieces

- |                |                                            |
|----------------|--------------------------------------------|
| <b>US-75KE</b> | transparent<br>extra thin: 0.05 mm         |
| <b>US-75KN</b> | transparent<br>standard thickness: 0.08 mm |



matrix band dispenser incl. matrix band roll, length: 15 m, width: 8.5 mm

- |                 |                                            |
|-----------------|--------------------------------------------|
| <b>US-15PKE</b> | transparent<br>extra thin: 0.05 mm         |
| <b>US-15PKN</b> | transparent<br>standard thickness: 0.08 mm |



matrix band refill roll, length: 15 m, width: 8.5 mm, package containing 3 rolls

- |                    |                                                                      |
|--------------------|----------------------------------------------------------------------|
| <b>US-15PKEN</b>   | transparent, extra thin: 0.05 mm                                     |
| <b>US-15PKNN</b>   | transparent, standard thickness: 0.08 mm                             |
| <b>US-15PKEN E</b> | package with 1 refill roll, transparent, extra thin: 0.05 mm         |
| <b>US-15PKNN E</b> | package with 1 refill roll, transparent, standard thickness: 0.08 mm |

