


preparation training

A-PTS

preparation trainer kit
mounting plate for frasco phantom heads with clip device for preparation trainers to exercise typical preparations
complete with: 2 square transparent plastic blocks printed with 13 different preparation designs in total, millimetric film for three-dimensional control and mounting key




A-PT

preparation trainer
2 square transparent plastic blocks printed with 13 different preparation designs in total, millimetric film for three-dimensional control



A-PTM

2 square plastic blocks simulating enamel, dentine and pulp layers made from 3 different layers of hard thermosetting frasco tooth material, printed with 13 different preparation designs in total, millimetric film for control




A-PT 99-001

2 square transparent plastic blocks printed with simple geometrical shapes "OSY" for preparation training, millimetric film for three-dimensional control

A-PTM 99-001

2 square plastic blocks simulating enamel, dentine and pulp layers printed with simple geometrical shapes „OSY“ for preparation training, millimetric film for control



A-PTSM

enamel, dentine, pulp preparation trainer kit mounting plate for frasco phantom heads with clip device for preparation trainers made from 3 different layers of hard thermosetting frasco tooth material
complete with: 2 square plastic blocks printed with 13 different preparation designs in total, millimetric film for control and mounting key



1531-SRM

mounting plate for screw-in or magnetic fixation of preparation or endodontic trainers in frasco phantom heads

1534-SRM

mounting plate used for fitting 2 frasco preparation trainers in one of the V-shaped positions left and right, intended to simulate the dental arch.
They enable optional screw-retention or retention using magnetic mounting plates in frasco phantom heads