

exodontic exercises

A-E M

extraction model for basic exercises
execution in cast aluminium, 32 metal teeth in CrNi-Steel,
embedded in resistant plastic resin, with back plate




A-E ZM

set of 32 metallic model teeth in CrNi-Steel for extraction
model A-E M

A-E ZME

individual metallic model teeth in CrNi-Steel
optional assortment



A-E K

extraction model for advanced exercises
execution in cast aluminium, 32 plastic teeth
embedded in resistant plastic resin, with back plate

plastic teeth roots will fracture in case of incorrect luxation




A-E ZK

set of 32 plastic model teeth in hard thermosetting material,
for extraction model A-E K

A-E ZKE

individual plastic model teeth in hard thermosetting material,
optional assortment



A-E MK

extraction model, execution in cast aluminium, 16 model teeth in
CrNi-Steel and 16 plastic teeth, embedded in resistant plastic resin,
with back plate

