

model teeth for pedodontal models

AK-6/2 ZPUL

Individual deciduous model tooth with screws, hollow pulp chamber for pulpotomy exercises, fitting dental models AK-6/2 and AK-6


upper	53, 54, 55
lower	85



AK-6/2 ZPUW

Individual deciduous model tooth with screws, pulp chamber filled with red wax for pulpotomy exercises e.g. with a spoon-excavator, fitting dental models AK-6/2 and AK-6


upper	53, 54, 55
lower	85



AK-6/2 ZPUWPK

Individual deciduous model tooth according to AK-6/2 ZPUW, pulp cavity filled with red wax, additional caries lesions suitable for models AK-6/2 and AK-6


upper	53, 54, 55
lower	85



AK-6/2 ZPUF

Individual deciduous model tooth for teeth endodontics suitable for model bases AK-6/2 and AK-6, screw retention, with transparent root and red hollow root canals, and ivory-coloured crowns made from Duroplast for pulpectomy exercises

upper	53, 54, 55
lower	85



AK-6/2 ZFK 16

Model tooth of the AK-6/2 series with fissure caries, upper 16

