

ANA-4 dental models and accessories

ANA-4 DAF	upper and lower model jaws combined with a demonstration articulator allowing simple functional movements (articulator fixed with screws to the model base).
ANA-4 V DAF	like ANA-4 DAF, however with simplified gingival contours
ANKA-4 DAF	upper and lower model jaws like ANA-4 DAF combined with a demonstration articulator, model teeth however not fixed with screws, but with the frasaco click system
ANKA-4 V DAF	like ANKA-4 DAF, however with simplified gingival countours
AN-4 WOK	elastic replacement gingiva for upper jaw, positioned over the model base, held in place by the model teeth, suitable for ANA-4 OK or ANKA-4 OK
AN-4 WOK KR	as above, however in red-transparent execution
AN-4 WOKV	as above, however in simplified execution suitable for ANA-4 V OK or ANKA-4 V OK
AN-4 WUK	elastic replacement gingiva for lower jaw, positioned over the model base, held in place by the model teeth, suitable for ANA-4 UK or ANKA-4 UK
AN-4 WUK KR	as above, however in red-transparent execution
AN-4 WUKV	as above, however in simplified execution suitable for ANA-4 V UK or ANKA-4 V UK
AN-4 WOKV 99-001	hard plastic replacement gingiva with simplified contours suitable for ANA-4 V OK or ANKA-4 V OK
AN-4 WUKV 99-001	hard plastic replacement gingiva with simplified contours suitable for ANA-4 V UK or ANKA-4 V UK
AN-4 RM SET	spare part: metal insert plate with thread (incl. screws) compatible with model series ANA-4 (V), ANKA-4 (V), AG-3, AG-3 KH, A-PB and A-PZ. It enables the optional screwing in or fastening with magnetic mounting plates in frasaco phantom heads.
AN-4 ROK	spare part: back plate for upper jaws of series ANA-4 (V) and ANKA-4 (V)
AN-4 RUK	spare part: back plate for lower jaws of series ANA-4 (V) and ANKA-4 (V)

