

DENTATUS LAMINEER TIPS FOR PROFIN/EVA

FLEXIBLE

The LTA-Flexible Great at Preserving Natural Esthetics

The Advantages

- Flexible design preserves natural esthetics
- Longer tip allows for even more accessibility
- For use in all Profin/EVA instruments

The new LTA-flexible tips are twice as flexible and 2mm longer than the standard LTA range of tips. These new features prevent ditching and iatrogenic damage to natural tooth and restorative contours in the anterior region raising the gold-standard in esthetic procedures.

Lamineer Tips can rotate freely to follow natural contours or be fixed for controlled detailing and shaping. When used in place of finishing strips in a free floating mode, the tips provide delicate control to shape and polish both the curved and flat surfaces to optimize esthetic form in the gingival, interproximal, facial and incisal areas.

More flexible and 2 mm longer

Premium diamond coating

Color coded



Dentatus LTA-Flexible, Lamineer Tips

Different Areas of Use



Dentatus AB
www.dentatus.com
SE: +46 8 546 509 00
info@dentatus.com
US: +1 212 481 1010
dentatus@dentatus.com

