DENTATUS IPR TIPS FOR PROFIN

FOR ALIGN ORTHODONTICS

Your Complete IPR Solution


The Advantages:

- Efficient and safe treatment
- Increased patient comfort
- Gentle to the gingiva
- One-hand operation

The solution consists of Interproximal Reduction (IPR) tips together with the reciprocating Profin PDX handpiece.



Profin PDX with IPR Tip



IPR in Align Orthodontics




Clear Aligner Treatment

Calibrated blade thickness from 0.1 mm to 0.5 mm

Long lasting tips

Autoclavable



Art.nr. LTA-P/SK2

Dentatus

Control and Performance Excellence with Profin's unique features

The Profin PDX handpiece has **unique features** such as reciprocating with a 1.2 mm axial stroke and self-steering tips. Alternatively the head has the possibility to lock the reciprocating tips in 8 different positions. Profin IPR's controlled reduction avoids rounding teeth thus preserving the natural anatomy of the tooth.IPR tips are color-coded to ensure time and performance excellence.

Product Description

You may choose to buy an IPR Starter Kit (**Art.nr. LTA-P/SK2**), including 5 tips and a stand or your chosen tips separately in a IPR refill pack of 2 pcs.

IPR Starter Kit



LTA-P/SK2



LTA-P3/2

IPR Refills Kits

All IPR Tips are color coded and have the grit size equivalent to the description below. Their color-coding is strictly related to blade thickness.

LTA-P1: Black — Double sided. Blade thickness 0.1 mm. Grit size 15µ.

LTA-P2: White — Double sided. Blade thickness 0.2 mm. Grit size 40µ.

LTA-P3: Yellow — Double sided. Blade thickness 0.3 mm. Grit size 50µ.

LTA-P4: Green — Double sided. Blade thickness 0.4 mm. Grit size 50µ.

LTA-P5: Red — Double sided. Blade thickness 0.5 mm. Grit size 50µ.

Accessory Tips

LTA-P25R: Violet — Safe-sided – abrasive on the right side. Blade thickness 0.25 mm. Grit size 50µ.

LTA-P25L: Violet — Safe-sided – abrasive on the left side. Blade thickness 0.25 mm. Grit size 50µ.

Safety Excellence

As safety is the ultimate concern for the dentist, we have developed two variants of our IPR tips (see product description below for details). IPR tips are either coated double-sided with diamond or safe-sided to ensure safety to adjacent teeth.

Dentatus AB
www.dentatus.com
SE: +46 8 546 509 00
info@dentatus.se
US: +1 212 481 1010
dentatus@dentatus.com

MADE IN SWEDEN

BL-6711-EN-04_2021-06