
Dentatus PROFIN® IPR LAMINEER® TIPS

PROFIN® IPR – FOR SAFE INTERPROXIMAL REDUCTION

At the Great New York Meeting in Chicago, November 2014, Dentatus pre-launched its newest addition in the range of Dentatus Lamineer tips for Profin/EVA. Tips for the highly acclaimed reciprocating Profin Directional System — PROFIN IPR™, for interproximal reduction.

PROFIN is an essential addition to any discipline in dentistry and, now more than ever, in orthodontics. For a dentist, patient safety is top priority; that is why we developed the PROFIN IPR. With PROFIN IPR, dentists can correct disproportionate tooth sizes, eliminate crowding, protect enamel and increase the durability of orthodontic treatment results. Profin IPR diminishes the risk of excessive stripping thereby ensuring more conservative removal of enamel in prescribed amounts.



THE ADVANTAGES:

- Safety is top concern for the patients
- Control is in your hand
- Performance made easier
- Esthetics is all natural

Safety is Top Concern for Our Patients
Profin IPR diminishes the risk of excessive stripping thereby ensuring more conservative removal of enamel in prescribed amounts. Profin IPR also emliminates the risk of damage to gingival tissue.

Control is in Your Hands
Used with the reciprocating Profin Directional System, Profin IPR´s movement is fluid as the 1.2 mm axial stroke allows for concise corrections and safe, quick and efficient modification


Performance Made Easier
Profin IPR´s controlled reduction avoids rounding teeth. It allows you to reach all the right places for optimal results, shaping, finishing and preserving the natural anatomy of the tooth.

Esthetics Have Never Been More Natural
With Profin IPR pacticioners can easily recreate teeth´s natural contours with a variety of Lamineer Tips in numerous shapes, sizes and grits.

Dentatus
MADE IN SWEDEN

Assortment:

Dentatus Profin LTA-P/SK, contains one of each tip below (number without /2) + a stand



Refills, 2 pcs/pkg:

LTA-P25R/2: violet shank, safe-sided – abrasive on right side when shank points down and back of tip towards user. Blade thickness 0.25mm

LTA-P25L/2: violet shank, safe-sided – abrasive on left side when shank points down and back of tip towards user. Blade thickness 0.25mm

LTA-P3/2: yellow shank. Double sided. Blade thickness 0.3mm

LTA-P4/2: green shank. Double sided. Blade thickness 0.4mm

LTA-P5/2: red shank. Double sided. Blade thickness 0.5mm

All IPR tips have the same grit – equivalent to 50µ. Color coding is relating to blade thickness only. No separate material safety data sheets available at this stage.

BL-6683-00_2015-02

Dentatus AB
Sweden +46 8 546 509 00, info@dentatus.se
www.dentatus.com
USA + 1 212 481 1010, dentatus@dentatus.com
MADE IN SWEDEN